

rapideye

Smart View Glasses
for intraoral scanner

Digital Solutions



Perfect
solution for
eye-hand
coordination


Easy
connection with
all intraoral
scanners

MODE MEDIKAL®

0-4-1-0-1-2

rapideye MG-40


Smart View Glasses
for intraoral scanner



High-Quality
Display

Ergonomic
Design

Simple
Connectivity



What is Smart View Glasses ?

Smart View Glasses offer a fantastic solution to minimize issues such as data loss and excessive data generation that may arise from the inability to maintain eye-hand coordination while using intraoral scanners. With these glasses, you can perform scans much more accurately and comfortably, thanks to the high-resolution image projected onto the glasses without the need to look at a computer screen.

Features

- Provides a full-screen detailed visual with transparency.
- Ergonomic design promotes correct work posture.
- Enables real-time synchronization with CAD systems.
- With specially designed nose pads, MG-40 Smart View Glasses can be comfortably worn over existing eyeglasses, providing a seamless and convenient experience for users with prescription eyewear.

System Requirements

- Wi-Fi Connection
- Windows 10 - Android 8 or higher

Works with all intraoral scanners



MODE MEDİKAL®

Abdi İpekçi Cad. No:58 Bayrampaşa
34030 İSTANBUL / TÜRKİYE
+90 (212) 612 64 09

